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| Test            | Cutoff (ng/mL) |
|-----------------|----------------|
| Amphetamine     | 500            |
| Methamphetamine | 500            |
| Morphine        | 300            |

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**Methamphetamine:** D-Methamphetamine (1000ng/mL), p-Hydroxy Methamphetamine (30,000ng/mL), L-Methamphetamine (8000ng/mL), Mephentermine (50,000ng/mL).

**Morphine:** Morphine (300ng/mL), Codeine (300ng/mL), Ethylmorphine (6250ng/mL), Hydrocodone (50,000ng/mL), Hydromorphone (3125ng/mL), Levorphanol (1500ng/mL), 6-Monoacetylmorphine (400ng/mL), Morphine 3-β-D-glucuronide (1000ng/mL), Oxycodone, (30,000ng/mL), Oxymorphone (100,000ng/mL).

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Acetaminophen, Acetophenetidin, N-Acetylprocainamide, Acetylsalicylic acid, Aminopyrine, Amoxicillin, Ampicillin, L-Ascorbic acid, Apomorphine, Aspartame, Atropine, Benzoic acid, Benzoic acid, Benzphetamine, Bilirubin, D/L-Brompheniramine, Caffeine, Cannabidiol, Chlorothiazide, Chloralhydrate, Chloramphenicol, D/L-Chloropheniramine, Chlorpromazine, Cholesterol, Chloroquine, Clonidine, Cortisone, L-Cotinine, Creatinine, Digoxin, β-Estradiol Diolofenac, Deoxycorticosterone, Dextromethorphan, Diflunisal, Diphenhydramine, Ecgonine methyl ester, Estrone-3-sulfate, L(-)-Epinephrine, Furosemide, [1R,2S] (-) Ephedrine, Ethyl-p-aminobenzoate, Erythromycin, Gentisic acid, Hemoglobin, Hydralazine, Hydrochlorothiazide, Hydrocortisone, O-Hydroxyhippuric acid, p-Hydroxyamphetamine, Iproniazid, p-Hydroxytyramine, D/L-Isoproterenol, Ibuprofen, Isoxsuprine, Ketamine, Ketoprofen, Labetalol, Loperamide, Meperidine, Meprobamate, Naloxone, Methoxyphenamine, Methylphenidate, Nalidixic acid, Naltrexone, Naproxen, Nicinamide, Nifedipine, Norethindrone, D-Norpropoxyphene, Noscaphine, D/L-Octopamine, D/L-Octopamine, Oxalic acid, Oxolinic acid, Oxymetazoline, Papaverine, Penicillin-G, Pentazocine hydrochloride, Perphenazine, Phenelzine, L-Phenylephrine, β-Phenylethylamine, Prednisolone, Prednisone, D/L-Propranolol, D-Propoxyphene, D-Pseudoephedrine, Quinacrine, Quinine, Quindine, Ranitidine, Salicylic acid, Serotonin, Sulfamethazine, Sulindac, Tetracycline, Thiamine, Thioridazine, D/L-Tyrosine, Tolbutamide, Triamterene, Trifluoperazine, Trimethoprim, Tryptamine, D/L-Tryptophan, Tyramine, Uric acid, Verapamil, Zomepirac, Tetrahydrocortisone 3 (β-D-glucuronide)

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- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 6th Ed. Biomedical Publ., Davis, CA. 2002; 129.
- FDA Guidance Document: Guidance for Premarket Submission for Kits for Screening Drugs of Abuse to be Used by the Consumer, 1997.
- A Handbook of Drug and Alcohol Abuse, Gail Winger, Third Edition, Oxford Press, 1992, page 146.
- Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735.
- Halows RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

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|   |                               |   |                         |   |                                    |
|---|-------------------------------|---|-------------------------|---|------------------------------------|
|  | Use by date                   |  | Caution                 |  | Catalog Number                     |
|  | Contains sufficient for tests |  | Manufacturer            |  | Do not reuse                       |
|  | Date of Manufacture           |  | Temperature limit       |  | Do not use if package is damaged   |
|  | Keep dry                      |  | Keep away from sunlight |  | In Vitro diagnostic medical device |
|  | Consult instructions for use  |  | Batch Code              |   |                                    |

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