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Agreement With Commercial Kit (%)		Agreement With GC/MS (%)	
Positive Agreement	98.8	Positive Agreement	98.6
Negative Agreement	100	Negative Agreement	97.5
Total Results	99.3	Total Results	98.0

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Drug Concentration Cut-Off Range	Rep	MOP	
		-	+
0% Cut-Off	30	30	0
-50% Cut-Off	30	30	0
-25% Cut-Off	30	28	2
Cut-Off	30	15	15
+25% Cut-Off	30	1	29
+50% Cut-Off	30	0	30

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Morphine (300ng/mL), Codeine (300ng/mL), Ethylmorphine (6250ng/mL), Hydrocodone (50000ng/mL), Hydromorphone (3125ng/mL), Levorphanol (1500ng/mL), 6-Monoacetylmorphine (400ng/mL), Morphine 3-β-D-glucuronide (1000ng/mL), Norcodeine (6250ng/mL), Normorphine (100000ng/mL), Oxycodone (30000ng/mL), Oxymorphone (100000ng/mL), Procaine (15000ng/mL), Thebaine (6250ng/mL),

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Acetaminophen, Acetophenetidin, N-Acetylprocainamide, Acetylsalicylic acid, Aminopyrine, Amoxicillin, Ampicillin, L-Ascorbic acid, Apomorphine, Aspartame, Atropine, Benzoic acid, Benzoic acid, Benzphetamine, Bilirubin, D/L-Brompheniramine, Caffeine, Cannabidiol, Chlorothiazide, Chloralhydrate, Chloramphenicol, D/L-Chlorpheniramine, Chlorpromazine, Cholesterol, Chloroquine, Clonidine, Cortisone, L-Cotinine, Creatinine, Digoxin, β-Estradiol Diclofenac, Deoxycorticosterone, Dextromethorphan, Diflunisal, Diphenhydramine, Egonine methyl ester, Estrone-3-sulfate, L(-)-Epinephrine, Furosemide, [1R,2S] (-) Ephedrine, Ethyl-p-aminobenzoate, Erythromycin, Gentic acid, Hemoglobin, Hydralazine, Hydrochlorothiazide, Hydrocortisone, O-Hydroxyhippuric acid, p-Hydroxyamphetamine, Iproniazid, p-Hydroxytyramine, D/L-Isoproterenol, Ibuprofen, Isoxsuprine, Ketamine, Ketoprofen, Labetalol, Loperamide, Meperidine, Meprobamate, Naloxone, Methoxyphenamine, Methylphenidate, Nalidixic acid, Naltrexone, Naproxen, Niacinamide, Nifedipine, Norethindrone, D-Norpropoxyphene, Noscapine, D/L-Octopamine, D/L-Octopamine, Oxalic acid, Oxolinic acid, Oxymetazoline, Papaverine, Penicillin-G, Pentazocine hydrochloride, Perphenazine, Phenelzine, L-Phenylephrine, β-Phenylethylamine, Prednisolone, Prednisone, D/L-Propranolol, D-Propoxyphene, D-Pseudoephedrine, Quinacrine, Quinine, Quindine, Ranitidine, Salicylic acid, Serotonin, Sulfamethazine, Sulindac, Tetracycline, Thiamine, Thioridazine, D/L-Tyrosine, Tolbutamide, Triamterene, Trifluoperazine, Trimethoprim, Tryptamine, D/L-Tryptophan, Tyramine, Uric acid, Verapamil, Zomepirac, Tetrahydrocortisone 3 (β-D-glucuronide)

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1. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 6th Ed. Biomedical Publ., Davis, CA. 2002; 129.
2. FDA Guidance Document: Guidance for Premarket Submission for Kits for Screening Drugs of Abuse to be Used by the Consumer, 1997.
3. A Handbook of Drug and Alcohol Abuse, Gail Winger, Third Edition, Oxford Press, 1992, page 146.
4. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735.
5. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

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	Use by date		Caution		Catalog Number
	Contains sufficient for tests		Manufacturer		Do not reuse
	Date of Manufacture		Temperature limit		Do not use if package is damaged
	Keep dry		Keep away from sunlight		In Vitro diagnostic medical device
	Consult instructions for use		Batch Code		

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